

Global Markets



Nicosh Automated External Defibrillator NC-05AED Plus

Compliance with ERC/AHA guidelines



NICOSH

Tel: +62 21-59987134 E-mail: ptnicoshindomedical@gmail.com
Fax: +62 21-59987134 website: nicoshmedlab.com
Kawasan Industri Milenium, Jalan Milenium Raya Blok K1 Nomor 1
Tangerang, Banten 15710 Indonesia

nicoshmedlab.com



NICOSH

Semi Automated External Defibrillators

can save life when a heart arrest happens

NC-05AED



Real-time battery status



Voice and visual instructions guide the rescuer through the entire defibrillation process



Smart Biphasic Truncated Exponential shock waveform



Remote management



7" HD screen (i5 only)



Portable & lightweight (2.3kg)



AED WALL CABINET



1 Turn ON the AED



2 Apply AED pads



3 Follow the voice instructions and deliver a shock



Features

- High capacity and long-life lithium manganese dioxide battery minimum 200 shocks 5-year life
- Visualized status of battery and system
- Detailed voice instructions and illustrations make rescue easier
- Automated tracking and direct data transfer via built-in 4G module (optional)
- Data storage of 24-hour ECG and up to 1000 events
- Deliver a 200J shock to adults and 50J shock to pediatric
- Daily self-tests report device status in time ideal for remote management